### IŞIKLA POLİMERİZE TEK KOMPONENTLİ DENTAL ADEZİV

#### 1. ÜRÜN BİLGİSİ

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### LIGHT POLYMERIZED SINGLE COMPONENT DENTAL ADHESIVE

#### 1. PRODUCT INFORMATION

For professional dentistry use only. It is a one-component, one-stage and one-application universal bonding agent that can be used in selection and total etching application, light cured with its excellent adhesion and adhesion strength, polymerized by radiation. The bonding agent is a nano-fillers and the film thickness is a 8 microns.

#### 2. PACKING

- TRIOBOND 5 ml in an easy-to-use unique bottle.

#### 3. CONTENTS

- 2-ethyl methacrylate, acetone, 4-morpholine, methacrylate monomer, silane, silica, water, initiator.

#### 4. USE

20 seconds for halogen/LED (700 mW/cm<sup>2</sup>), high energy LED (1200 mW/cm<sup>2</sup>) 5 – 10 seconds. If you are storing TRIOBOND in the refrigerator, let the room temperature for a few minutes.

#### 5. INDICATIONS

- Direct bonding of light-cured composites and compomers to tooth structure,
- Repair of composite, porcelain zirconium etc. restorations inside and outside the mouth,
- Treatment of hypersensitive teeth,
- Sealing of tooth preparation for indirect restorations.

#### 6. CONTRAINDICATIONS

- Pulp capping,
- Avoid use of this product in patients with known allergies,
- With chemically curing composites

#### 7. WARNINGS

- When a problem is encountered in allergic and hypersensitive persons, the application should be stopped, washed with plenty of water and a specialist doctor should be consulted.
- Personal protection such as gloves, face mask and safety glasses should be worn at all time.
- Avoid inhalation or ingestion of material.
- TRIOBOND is flammable and volatile. Do not use directly next to flammable. Keep away from explosive sources. Keep away from direct sunlight.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Do not mix with products such as eugenol, diamine silver fluoride, parachlorophenol, guaiacol, phenols, aluminum chloride, ferric sulphate, aluminum sulphate.
- Dispose of contents and bottle according to local regulations.

#### 8. INSTRUCTIONS FOR APPLICATION

- Prepare the tooth using standard techniques. Please be sure to thoroughly clean the temporary material from the surface. Use rubber-dam for protection. Wash and dry the prepared tooth surfaces by lightly applying air with an air syringe.
- TRIOBOND is a self-etch bond with superior bonding strength and is suitable for use with total-etch and selective etch techniques.
- Before use, shake the bottle thoroughly, open the cap of the bottle, apply a few drops of it to the bowl by applying gentle pressure from the notches on the sides of the bottle.
- The TRIOBOND with the applicator should be applied to the tooth as soon as possible, while waiting may cause adverse effects in terms of adhesion strength.
- After the TRIOBOND is taken into the body, the cap of the bottle should be closed immediately.
- Wait 10 seconds after applying the TRIOBOND to the cavity, Dry the cavity with air spray for 5 seconds so that the excess TRIOBOND is removed from the cavity. The air coming from the compressor must be dry. The moist and oily incoming air negatively affects the adhesion


strength of the TRIOBOND. When using on vital teeth, do not dry the applied area completely. Residue etc. on the tooth surface. If there are substances, the surface should be cleaned as follows. It is cleaned with alcohol or citric acid on the tooth surface.

- Irradiating for halogen/LED (700 mW/cm<sup>2</sup>) for 20 seconds, high energy LED (1200 mW/cm<sup>2</sup>) 5 – 10 seconds provides sufficient polymerization.
- Keep the light source as close as possible.
- Protect from the harmful effects of light with protective glasses.
- Bond removed form the bottle should be quickly.
- If the applied material has to come into contact with water, blood or saliva before it polymerizes, wash it, dry it and repeat the procedure by applying the material again.
- For the treatment of hypersensitive teeth, the self-etch technique is recommended.

#### 9. STORAGE CONDITIONS

Store the bottle at 1 to 25 centigrade degrees. The shelf life of the product is 2 years. For repeated use, make sure that the cap is closed correctly after use. After use, the product should be stored in a cool, dry and dark environment.

#### 10. SYMBOLS

|   |   |
|---|---|
|    | Date of production  |
|    | Manufacture   |
|    | Refer to instruction of use                                       |
|    | Lot number  |
|    | Temperature limit   |
|   | Caution   |
|  | CE Declaration with Notified Body Number                          |
|  | Keep away from sunlight.  |
|  | Keep dry.   |
|  | Product reference code  |
|  | Medical Device  |
|  | Professional for use, only.                                       |
|  | Flammable.  |
|  | Do not use if package is damaged ad consult instructions for use. |